

CD301

Socket Extractor



Customer Service
(800) 819-5980 or (208) 429-0026
www.coyotedesign.com



Advena Ltd
Pure Offices, Plato Close
Tachbrook Park
Warwick, CV34 6WE, UK

MADE IN USA



GD301.rev.07152014

The threaded rod and wheel both come pre-lubricated, however it is a good idea to lubricate these items about every other month.

Before pouring mold drill a 1/4" hole into the proximal end of the pipe you are going to use in the mold.

With bench vice, clamp the bottom of the extractor frame.



**A. Extractor Frame Top B. U-Bolt C. 1/4 Hex Head Bolts (2)
 D. 1/4 Nuts E. U-Bolt Attachment Screws (4) F. Frame Bolt
 (4) G. Coupler H. Hand Wheel I. Frame & Coupler Washers
 (12) J. Extractor Rod K. Threaded Rod L. Frame Cap M.
 Left Frame Arm N. Right Frame Arm**


Correctly Using the Coyote Design Socket Extractor



**Slide right and left frame arms into top extractor frame.
 Using provided nuts, bolts, and washers, screw both arms into place on top frame. Ensure washers are on both sides of frame.
 Using allen wrench and open wrench to tighten nuts.
 Drop threaded rod into bushing.**



**Thread hand wheel on to threaded rod.
 Screw in handles.
 Attach coupler to threaded rod with supplied screws.
 Slide extractor rod through frame arm holes and through extractor U bolt.**



**Trim as close to the trim lines on the socket as possible.
 Slide coupler over threaded rod and insert provided bolt into hole. Insert second bolt into socket pipe hole.
 Using provided screws attach the anchor to the distal end of the socket to be pulled.
 Slide anchor rod through extractor frame and anchor. Turn wheel clockwise to extract mold.
*Retro fitable for pediatric sockets.***